

```

1 # Weeek4 - Homework : Chest X-Ray Image Classification (Normal or Pneumonia) using Transfer Learning
2 #1- Download Chest X-Ray Image Data From Kaggle
3 #2- Plot Normal and Pneumonia Chest X-Ray Images
4 #3- Chest X-Ray Image Preprocessing/ Augmentation/ Transformation for Model Training, Validation and Testing Dataset
5 #4- Create InceptionV3 Base CNN Model
6 #5- Create a Sequential Model for traing & prediction. Add base model, pooling layer & fully connected Dense layers
7 #6- Compile the Sequential Model and Display Model Summary
8 #7- Fit and Train the Sequential Model and Display Results
9 #8- Evaluate the Sequential Model Performance and Show Results
10 #9- Test the pre-trained Sequential Model and Show Results
11 # Work below is based on Anjana Tiha github project at https://github.com/anjanatiha/Pneumonia-Detection-from-Chest-X-Ray-Images-with-D
12 import os #for OS dependent functionalities
13 from os import listdir, makedirs, getcwd, remove
14 import numpy as np # data presprocessing
15 import pandas as pd
16 import matplotlib
17 from matplotlib import pyplot as plt # visualize the data
18 import matplotlib.image as ming
19 import imgaug.augmenters as iaa
20 import math, json
21 import datetime
22 import random
23 import re
24 import cv2
25 from PIL import Image #Python Imaging Library
26 from pathlib import Path
27 from glob import glob # find pathnames for images matching specified pattern
28 import keras
29 # Inception-v3 is convolution neural network trained on million+ images from ImageNet db.
30 # The network is 48 layers deep and can classify images into 1000 object categories
31 # The Inception-v3 model expects color images to have the square shape 299x299.
32 from keras.applications.inception_v3 import InceptionV3
33 from keras.preprocessing.image import ImageDataGenerator #to perform image augmentation on the fly & easy way
34 from keras.models import Sequential
35 from keras.layers import Conv2D, MaxPooling2D
36 from keras.layers import Activation, Dropout, Flatten, Dense
37 from keras.utils import plot_model
38 from keras import optimizers
39 from keras.optimizers import Adam, SGD , RMSprop
40 from keras import backend as K
41 import warnings
42 warnings.filterwarnings('ignore')
43 import tensorflow as tf
44 import torch # to save model etc on google colab drive
45 from google.colab import drive # to save model etc on google colab drive
46
47 #!pip install tensorflow --upgrade
48 #!pip install keras --upgrade
49
50 print("Tensorflow version " + tf.__version__)
51
52 print("Keras version " + keras.__version__)
53
54 print("Numpy version " + np.__version__)
55
56
57 ##### Download Chest X-Ray Image Data From Kaggle #####
58
59 print("\n" + '\033[1m \033[4mDownload Chest X-Ray Data From Kaggle\033[0m' + "\n")
60 # Setup the Kaggle environemt for image download
61 #os.environ['KAGGLE_USERNAME'] = 'XXXX'
62 #os.environ['KAGGLE_KEY'] = 'XXXXXXX'
63 # Download the Kaggle Chest X-Ray dataset
64 #!kaggle datasets download -d paultimothymooney/chest-xray-pneumonia
65 # unzip the chest x-ray images
66 #!unzip chest-xray-pneumonia.zip
67 #!unzip chest_xray.zip
68
69 ##### Plot Normal and Pneumonia Chest X-Ray Images #####
70
71 print("\n" + '\033[1m \033[4mPlot Normal and Pneumonia Chest X-Ray Images\033[0m' + "\n")
72
73 # Define and configure file directories for train, validate and test image dataset
74 training_dir = './chest_xray/train'
75 validate_dir = './chest_xray/val'
76 test_dir = './chest_xray/test'
77
78 # Plot an image of Normal Chest Xray
79 plt.figure(1, figsize = (12 , 6))
80 plt.subplot(1 , 2 , 1)
81 print('Number of Normal Images: ' + str(len(glob(training_dir+"/NORMAL/*.jpeg"))))
82 img = np.asarray(plt.imread(glob(training_dir+"/NORMAL/*.jpeg")[0]))
83 plt.title('Normal Chest X-Ray')
84 plt.imshow(img)
85
86 # Plot an image of Pneumonia Chest Xray
87 plt.subplot(1 , 2 , 2)
88 img = np.asarray(plt.imread(glob(training_dir+"/PNEUMONIA/*.jpeg")[0]))

```

```

89 print('Number of Pneumonia Images: ' + str(len(glob(training_dir+"/PNEUMONIA/*.jpeg"))))
90 print("")
91 plt.title('Pneumonia Chest X-Ray')
92 plt.imshow(img)
93 plt.show()
94
95 ##### Chest X-Ray Image Preprocessing/ Augmentation/ Transformation for Model Training, Validation and Testing Dataset #####
96
97 print("\n" + '\033[1m \033[4mImage Preprocessing/ Augmentation/ Transformation for Training, Validation and Testing Dataset\033[0m' + "
98
99 # Define Constants used
100 # Rescale the pixel values (between 0 and 255) to [0, 1] interval
101 # Neural networks perform better with normalize data.
102 # rescale normalizes the image pixel values to have zero mean & standard deviation of 1.
103 rescale = 1./255 # Scale the image between 0 and 1
104 target_size = (150, 150) # by keras doc, conv2d input shape should be (height, width) == (row, column) sequence
105 batch_size=32
106 img_height=150
107 img_width=150
108
109 # Use ImageDataGenerator() so Python generators that automatically turn image files into preprocessed tensors that can be fed directly
110 #1. Decode the JPEG content to RGB grids of pixels.
111 #2. Convert these into floating-point tensors.
112 #3. Rescale the pixel values (between 0 and 255) to the [0, 1] interval ( neural networks perform better with normalize data).
113 #4. Helps us easily augment images.
114
115 # Create Training Image Data Generator and augment the training images to augment our data-set and improve generalization.
116 training_datagen = ImageDataGenerator(
117     rescale=rescale,
118     width_shift_range=0.2,
119     height_shift_range=0.2,
120     shear_range=0.2,
121     zoom_range=0.2,
122     horizontal_flip=True
123 )
124
125 print('Training Data Generation')
126 # this is the augmentation configuration we will use for training
127 training_generator = training_datagen.flow_from_directory(
128     training_dir,
129     target_size=target_size,
130     class_mode='binary',
131     batch_size=32
132 )
133 print('Training classes are' + str(training_generator.class_indices) + '\n')
134
135
136 # Create Vaildation Image Data Generator
137 # No image data augmentation for validation image data generation, only rescale
138 validation_datagen = ImageDataGenerator(rescale=rescale)
139 print('\n' + 'Validation Data Generation')
140 # this is the augmentation configuration we will use for validation
141 validation_generator = validation_datagen.flow_from_directory(
142     validate_dir,
143     target_size=target_size,
144     class_mode='binary',
145     batch_size=32
146 )
147 print('Validation classes are' + str(validation_generator.class_indices) + '\n')
148
149 # Create Testing Image Data Generator
150 test_datagen = ImageDataGenerator(rescale=rescale)
151 print('\n' + 'Test Data Generation')
152 # this is the augmentation configuration we will use for testing
153 test_generator = test_datagen.flow_from_directory(
154     test_dir,
155     target_size=target_size,
156     class_mode='binary',
157     batch_size=1
158 )
159 print('Test classes are' + str(test_generator.class_indices))
160
161 ##### Create InceptionV3 Base CNN Model #####
162
163 print("\n" + '\033[1m \033[4mCreating Pre-Trained InceptionV3 Base Model For Transfer Learning\033[0m' + "\n")
164
165 # A Convolutional Neural Network (CNN) architecture has three main parts:
166 # A convolutional layer that extracts features from a source image. Convolution helps with blurring, sharpening, edge detection, noise
167 # A pooling layer that reduces the image dimensionality without losing important features or patterns.
168 # A fully connected layer also known as the dense layer, in which the results of the convolutional layers are fed through one or more n
169
170 # Create the base pre-trained InceptionV3 CNN model for image processing
171 # "include_top" argument set to False so fully-connected output layers of the model
172 # used to make predictions is not loaded, allowing a new output layer to be added and trained.
173 # A model without a top will output activations from the last convolutional or pooling layer directly.
174 # Your input should be input_shape=(img_rows, img_cols, 3 channels RGB) or (batchSize, height, width, feature/channels)
175 base_model = InceptionV3(weights='imagenet', include_top=False, input_shape=(150, 150, 3))
176
177 ##### Create a Sequential Model for traing & prediction. Add base model, pooling layer & fully connected Dense layer
178 # Initializing the network using the Sequential Class

```

```

179 # The Sequential model (a linear stack of layers) is a prediction model which is trained with a set of training sequences.
180 # Once trained, the sequential model is used to perform sequence predictions.
181 # Sequential class which is a linear stack of Layers where after creating you can define all of the layers in the constructor
182 model = Sequential()
183 # Add the convolution layer by adding the pre-trained InceptionV3 base model output
184 model.add(base_model)
185 # Add a pooling layer by adding global spatial average pooling layer
186 model.add(keras.layers.GlobalAveragePooling2D())
187 # Add a fully-connected layer;
188 model.add(Dense(1024, activation='relu')) # rectified linear unit
189 # Dropout consists in randomly setting a fraction rate of input units to 0 at each update during training time, which helps prevent over
190 model.add(Dropout(0.5))
191 # Normalize the activations of the previous layer at each batch, i.e. applies a transformation that maintains the mean activation close
192 model.add(keras.layers.BatchNormalization())
193 # Model output arrays will be of shape (*, 8)
194 model.add(Dense(8, activation='softmax'))
195 # Add logistic layer
196 # Using 1 class for binary classification with activation sigmoid
197 # Model output arrays will be of shape (*, 1)
198 model.add(Dense(1, activation='sigmoid'))
199
200 ##### Compile the Model and Display Model Summary #####
201
202 print('Compiling The Model' + "\n")
203
204 model.compile(loss= 'binary_crossentropy',
205               optimizer='rmsprop',          # root mean square propagation
206               metrics=['accuracy'])
207
208 print('Displaying The Model Summary Information below' + "\n")
209 model.summary()
210 print("")
211
212 ##### Fit and Train the Sequential Model and Display Results #####
213
214 print("\n" + '\033[1m \033[4mStarting to Fit and Train The Model\033[0m' + "\n")
215
216 # Train the model on data generated batch-by-batch by a Python generator
217 # Train the model for defined epochs for steps_per_epoch times in each epoch
218 # An epoch usually means one iteration over all of the training data.
219 # For instance for 20,000 images and a batch size of 100, the epoch should contain 20,000 / 100 = 200 steps.
220
221 steps_per_epoch=len(training_generator)
222 print('steps_per_epoch is ' + str(steps_per_epoch) + '\n')
223 validation_steps=len(validation_generator)
224 print('validation_steps is ' + str(validation_steps) + '\n')
225
226 history = model.fit_generator(
227     training_generator,
228     steps_per_epoch=steps_per_epoch,
229     epochs=5,
230     verbose=1,
231     validation_data=validation_generator,
232     validation_steps=validation_steps,
233 )
234
235 # Save the Weights and the Model
236 print("")
237 drive.mount('/content/gdrive') # if drive is dismounted then will need to forget previous authentication key and reauthenticate
238 model_save_name = 'ChestXRayClassifier.pt'
239 model_path_filename = F"/content/gdrive/My Drive/{model_save_name}"
240 print('\n' + 'Saving the trained model at %s ' % path)
241 torch.save(model, model_path_filename) #torch.save(model.state_dict(), path) gives error
242
243 # Plot the Model and its layers
244 path = F"/content/gdrive/My Drive"
245 plot_model(model, to_file='chest-xray-classification-model.png')
246 model_img = Image.open("chest-xray-classification-model.png")
247 plt.figure(1, figsize = (16 , 16))
248 plt.title('Trained Model')
249 plt.imshow(model_img)
250
251 print('\n' "Completed Model Fitting and Training" '\n')
252
253 ##### Evaluate Model Performance #####
254
255 print("\n" + '\033[1m \033[4mPerform and Plot Model Performance Evaluation\033[0m' + "\n")
256
257 xlabel = 'Epoch'
258 legends = ['Training', 'Validation']
259 ylim_pad = [0.01, 0.1]
260
261 plt.figure(figsize=(15, 5))
262
263 # Plot training & validation Accuracy values
264
265
266 y1 = history.history['acc']
267 y2 = history.history['val_acc']
268

```

```

269 min_y = min(min(y1), min(y2))-ylim_pad[0]
270 max_y = max(max(y1), max(y2))+ylim_pad[0]
271
272 plt.subplot(121)
273
274 plt.plot(y1)
275 plt.plot(y2)
276
277 plt.title('Model Accuracy', fontsize=17)
278 plt.xlabel(xlabel, fontsize=15)
279 plt.ylabel('Accuracy', fontsize=15)
280 plt.ylim(min_y, max_y)
281 plt.legend(legends, loc='upper left')
282 plt.grid()
283
284 # Plot training & validation loss values
285
286 y1 = history.history['loss']
287 y2 = history.history['val_loss']
288
289 min_y = min(min(y1), min(y2))-ylim_pad[1]
290 max_y = max(max(y1), max(y2))+ylim_pad[1]
291
292 plt.subplot(122)
293
294 plt.plot(y1)
295 plt.plot(y2)
296
297 plt.title('Model Loss', fontsize=17)
298 plt.xlabel(xlabel, fontsize=15)
299 plt.ylabel('Loss', fontsize=15)
300 plt.ylim(min_y, max_y)
301 plt.legend(legends, loc='upper left')
302 plt.grid()
303
304 plt.show()
305
306 ##### Test the pre-trained Sequential Model and Show Results #####
307
308 print("\n" + '\033[1m \033[4mTesting The Trained Model\033[0m' + "\n")
309
310 test_steps=len(test_generator)
311 results = model.evaluate_generator(generator=test_generator, steps=test_steps, verbose=1)
312 print('\n' + 'Trained Model Testing Results: Loss and Accuracy: ' + str(results))
313
314 history = model.fit_generator(
315     training_generator,
316     steps_per_epoch=steps_per_epoch,
317     epochs=5,
318     verbose=1,
319     validation_data=validation_generator,
320     validation_steps=validation_steps,
321 )
322
323 # Save the Weights and the Model
324 print("")
325 drive.mount('/content/gdrive') # if drive is dismounted then will need to forget previous authentication key and reauthenticate
326 model_save_name = 'ChestXRayClassifier.pt'
327 model_path_filename = F"/content/gdrive/My Drive/{model_save_name}"
328 print('\n' + 'Saving the trained model at %s ' % path)
329 torch.save(model, model_path_filename) #torch.save(model.state_dict(), path) gives error
330
331 # Plot the Model and its layers
332 path = F"/content/gdrive/My Drive"
333 plot_model(model, to_file='chest-xray-classification-model.png')
334 model_img = Image.open("chest-xray-classification-model.png")
335 plt.figure(1, figsize = (16 , 16))
336 plt.title('Trained Model')
337 plt.imshow(model_img)
338
339 print('\n' "Completed Model Fitting and Training" '\n')
340
341 ##### Evaluate Model Performance #####
342
343 print("\n" + '\033[1m \033[4mPerform and Plot Model Performance Evaluation\033[0m' + "\n")
344
345 xlabel = 'Epoch'
346 legends = ['Training', 'Validation']
347
348 ylim_pad = [0.01, 0.1]
349
350 plt.figure(figsize=(15, 5))
351
352 # Plot training & validation Accuracy values
353
354 y1 = history.history['acc']
355 y2 = history.history['val_acc']
356
357 min_y = min(min(y1), min(y2))-ylim_pad[0]
358 max_y = max(max(y1), max(y2))+ylim_pad[0]

```

```

359 plt.subplot(121)
360 plt.plot(y1)
361 plt.plot(y2)
362 plt.title('Model Accuracy', fontsize=17)
363 plt.xlabel(xlabel, fontsize=15)
364 plt.ylabel('Accuracy', fontsize=15)
365 plt.ylim(min_y, max_y)
366 plt.legend(legends, loc='upper left')
367 plt.grid()
368
369 # Plot training & validation loss values
370
371 y1 = history.history['loss']
372 y2 = history.history['val_loss']
373
374 min_y = min(min(y1), min(y2))-ylim_pad[1]
375 max_y = max(max(y1), max(y2))+ylim_pad[1]
376
377 plt.subplot(122)
378 plt.plot(y1)
379 plt.plot(y2)
380
381 plt.title('Model Loss', fontsize=17)
382 plt.xlabel(xlabel, fontsize=15)
383 plt.ylabel('Loss', fontsize=15)
384 plt.ylim(min_y, max_y)
385 plt.legend(legends, loc='upper left')
386 plt.grid()
387
388 plt.show()
389
390 ##### Test the trained Model #####
391
392 print("\n" + '\033[1m \033[4mTesting The Trained Model\033[0m' + "\n")
393
394 test_steps=len(test_generator)
395 results = model.evaluate_generator(generator=test_generator, steps=test_steps, verbose=1)
396 print('\n' + 'Trained Model Testing Results: Loss and Accuracy: ' + str(results)).
400

```



Tensorflow version 1.14.0
Keras version 2.2.5
Numpy version 1.16.5

Download Chest X-Ray Data From Kaggle

Plot Normal and Pneumonia Chest X-Ray Images

Number of Normal Images: 1341
Number of Pneumonia Images: 3875




Image Preprocessing/ Augmentation/ Transformation for Training, Validation and Testing Dataset

Training Data Generation
Found 5216 images belonging to 2 classes.
Training classes are{'NORMAL': 0, 'PNEUMONIA': 1}

Validation Data Generation
Found 16 images belonging to 2 classes.
Validation classes are{'NORMAL': 0, 'PNEUMONIA': 1}

Test Data Generation
Found 624 images belonging to 2 classes.
Test classes are{'NORMAL': 0, 'PNEUMONIA': 1}

Creating Pre-Trained InceptionV3 Base Model For Transfer Learning

Compiling The Model

Displaying The Model Summary Information below

Model: "sequential_17"

Layer (type)	Output Shape	Param #
=====		
inception_v3 (Model)	(None, 3, 3, 2048)	21802784
global_average_pooling2d_17	(None, 2048)	0
dense_49 (Dense)	(None, 1024)	2098176
dropout_17 (Dropout)	(None, 1024)	0
batch_normalization_1787 (Ba	(None, 1024)	4096
dense_50 (Dense)	(None, 8)	8200
dense_51 (Dense)	(None, 1)	9
=====		
Total params: 23,913,265		
Trainable params: 23,876,785		
Non-trainable params: 36,480		

Starting to Fit and Train The Model

steps_per_epoch is 163

validation_steps is 1

Epoch 1/5

163/163 [=====] - 836s 5s/step - loss: 0.5630 - acc: 0.7519 - val_loss: 0.8313 - val_acc: 0.5000

Epoch 2/5

163/163 [=====] - 760s 5s/step - loss: 0.5065 - acc: 0.7837 - val_loss: 0.9058 - val_acc: 0.5000

Epoch 3/5

163/163 [=====] - 756s 5s/step - loss: 0.4917 - acc: 0.7868 - val_loss: 0.7002 - val_acc: 0.5000

Epoch 4/5

163/163 [=====] - 759s 5s/step - loss: 0.4371 - acc: 0.8206 - val_loss: 0.5817 - val_acc: 0.6875

Epoch 5/5


163/163 [=====] - 765s 5s/step - loss: 0.4001 - acc: 0.8430 - val_loss: 0.5938 - val_acc: 0.6875



Drive already mounted at /content/gdrive; to attempt to forcibly remount, call drive.mount("/content/gdrive", force_remount=True)

Saving the trained model at /content/gdrive/My Drive

Completed Model Fitting and Training

Perform and Plot Model Performance Evaluation





Testing The Trained Model

624/624 [=====] - 60s 96ms/step

Trained Model Testing Results: Loss and Accuracy: [0.5722941191723714, 0.6634615384615384]

Epoch 1/5

26/163 [==>.....] - ETA: 11:26 - loss: 0.4076 - acc: 0.8269

```
KeyboardInterrupt                                Traceback (most recent call last)
<ipython-input-22-a73bc0487cf3> in <module>()
    307     verbose=1,
    308     validation_data=validation_generator,
--> 309     validation_steps=validation_steps,
    310 )
    311
```

⏏ 6 frames

```
/usr/local/lib/python3.6/dist-packages/tensorflow/python/client/session.py in __call__(self, *args, **kwargs)
    1456     ret = tf_session.TF_SessionRunCallable(self._session._session,
    1457                                           self._handle, args,
-> 1458                                           run_metadata_ptr)
    1459     if run_metadata:
    1460         proto_data = tf_session.TF_GetBuffer(run_metadata_ptr)
```

KeyboardInterrupt:

SEARCH STACK OVERFLOW

